



Product Datasheet

ANXA1 antibody (orb125350)

Description	Rabbit polyclonal antibody to ANXA1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04083
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 2...



Western blot analysis of extracts from J...