


Product Datasheet

53BP1 antibody (orb125356)

Description	Rabbit polyclonal antibody to 53BP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q12888
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Peptide - +

Immunohistochemical analysis of formalin...