



Product Datasheet

p47 phox antibody (orb125362)


Description	Rabbit polyclonal antibody to p47 phox
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P14598
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from C...



Western blot analysis of extracts from H...



Confocal immunofluorescence analysis of ...