



## Product Datasheet

### Cyclin E1 antibody (orb125365)


|                            |                                                                           |
|----------------------------|---------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Cyclin E1                                   |
| <b>Species/Host</b>        | Rabbit                                                                    |
| <b>Reactivity</b>          | Human                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                              |
| <b>Tested Applications</b> | IF, IHC, WB                                                               |
| <b>Immunogen</b>           | C-terminal                                                                |
| <b>Concentration</b>       | 1mg/ml                                                                    |
| <b>Note</b>                | For research use only                                                     |
| <b>Clonality</b>           | Polyclonal                                                                |
| <b>Uniprot ID</b>          | <b>P24864</b>                                                             |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:500-1:1000, ELISA: 1:10000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                           |



Western blot analysis of extracts from H...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...