

Product Datasheet

HDAC4 HDAC5 HDAC9 Antibody (orb1253696)

Description

HDAC4 HDAC5 HDAC9 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

HDAC4/HDAC5/HDAC9 (phospho-Ser246/259/220) antibody was raised against a peptide sequence around phosphorylation site of serine 246/259/220 (T-A-S (p) -EP) derived from Human HDAC4/HDAC5/HDAC9.

Target

HDAC4 HDAC5 HDAC9

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000,
Immunohistochemistry: 1:50~1:100

Clonality

Polyclonal

MW

140, 124 kDa

Uniprot ID

[P56524](#), [Q9UKV0](#), [Q9UQL6](#)

NCBI


[NP_055522.1](#), [NP_001015053.1](#),
[NP_006028.2](#)

Dilution Range


Western Blot: 1:500~1:1000,
Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.



Western blot analysis
of lysed extracts ...



Immunohistochemical
analysis of paraffin...