

Product Datasheet

RPS6KA1 Antibody (orb1253792)

Description

RPS6KA1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IF, IHC, WB

Applications

Immunogen

p90RSK (Phospho-Thr348) antibody was raised against a peptide sequence around phosphorylation site of threonine 348 (S-R-T (p) -P-R) derived from Human p90RSK.

Target

RPS6KA1

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000,
Immunohistochemistry: 1:50~1:100,
Immunofluorescence: 1:100~1:200

Clonality

Polyclonal

MW

90 kDa

Uniprot ID

[Q15418](#)

NCBI


[NP_001006666.1](#)

Dilution Range


Western Blot: 1:500~1:1000,
Immunohistochemistry: 1:50~1:100,
Immunofluorescence: 1:100~1:200

Expiration Date

12 months from date of receipt.



Western blot analysis
of lysed extracts ...



Top Image:
Immunohistochemical
analysis ...