



Product Datasheet

EGFR antibody (orb12538)


Description	Rabbit polyclonal antibody to EGFR
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human EGFR(NP_005219.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	175KDa
Uniprot ID	P00533
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC12 cell...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of HeLa cell...