


Product Datasheet

SLC12A2 Antibody (orb1253857)

Description

SLC12A2 Antibody

Species/Host	Gallus
Reactivity	Bovine, Canine, Gallus, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	800-1000
Target	SLC12A2
Preservatives	Phosphate-Buffered Saline. No preservatives added.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SLC12A2 antibody can be stored at 4°C for short term (weeks). Long term storage should be at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Solute carrier family 12 member 2 antibody can be used for the detection of Solute carrier family 12 member 2 by Western Blot.
Clonality	Polyclonal
MW	131.5 kDa (calculated)
Uniprot ID	P55011
NCBI	NP_001037.1
Dilution Range	Solute carrier family 12 member 2 antibody can be used for the detection of Solute carrier family 12 member 2 by Western Blot.
Expiration Date	12 months from date of receipt.



E coli-derived fusion protein as test an...