

Product Datasheet

SLC9A3R1 Antibody (orb1253989)

Description

SLC9A3R1 Antibody

Species/Host

Gallus

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Direct Elisa
Test. Protein
as test
antig...

Immunogen

240-358

Target

SLC9A3R1

Preservatives

Phosphate-Buffered Saline. No preservatives added.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SLC9A3R1 antibody can be stored at 4°C for short term (weeks). Long term storage should be at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Ezrin-radixin-moesin binding phosphoprotein 50 antibody can be used for the detection of Ezrin-radixin-moesin binding phosphoprotein 50 by ELISA and Western Blot, may also work for ICC and IHC.

Clonality

Polyclonal

MW

38.7 kDa (calculated)

Uniprot ID

[O14745](#)

NCBI

[NP_004243.1](#)

Dilution Range

Ezrin-radixin-moesin binding phosphoprotein 50 antibody can be used for the detection of Ezrin-radixin-moesin binding phosphoprotein 50 by ELISA and Western Blot, may also work for ICC and IHC.

Expiration Date

12 months from date of receipt.

