

Product Datasheet

IL8 Antibody (orb1254052)

Description

IL8 Antibody

Species/Host

Gallus

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

28-99

Target

IL8

Preservatives

Phosphate-Buffered Saline. No preservatives added.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

IL8 antibody can be stored at 4°C for short term (weeks). Long term storage should be at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Interleukin-8 antibody can be used for the detection of Interleukin-8 by ELISA and Western Blot, may also work for ICC and IHC.

Clonality

Polyclonal

MW

11.1 kDa (calculated)

Uniprot ID

[P10145](#)

NCBI


[P10145](#)

Dilution Range

Interleukin-8 antibody can be used for the detection of Interleukin-8 by ELISA and Western Blot, may also work for ICC and IHC.

Expiration Date

12 months from date of receipt.



Recombinant protein as test antigen (5 u...