



Product Datasheet

AKAP14 antibody (orb125406)

Description	Rabbit Polyclonal to AKAP14.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q86UN6
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from J...



Confocal immunofluorescence analysis of ...