

## Product Datasheet

### GC Antibody (orb1254065)

## Description

GC Antibody

### Species/Host

Gallus

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

### Immunogen

1-474

### Target

GC

### Preservatives

Phosphate-Buffered Saline. No preservatives added.

### Form/Appearance

Liquid

### Concentration

1 mg/mL

### Storage

GC antibody can be stored at 4°C for short term (weeks). Long term storage should be at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

Vitamin D-binding protein antibody can be used for the detection of Vitamin D-binding protein by ELISA and Western Blot, may also work for ICC and IHC.

### Clonality

Polyclonal

### MW

53 kDa (calculated)

### Uniprot ID

[P02774](#)

### NCBI


[NP\\_000574.2](#)

### Dilution Range

Vitamin D-binding protein antibody can be used for the detection of Vitamin D-binding protein by ELISA and Western Blot, may also work for ICC and IHC.

### Expiration Date

12 months from date of receipt.



Human plasma Gc-globulin protein as test...