


Product Datasheet

AKR1CL1 antibody (orb125410)

Description	Rabbit Polyclonal to AKR1CL1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q5T2L2
Dilution Range	IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...