


Product Datasheet

ALPK2 antibody (orb125415)

Description	Rabbit Polyclonal to ALPK2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q86TB3
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...