



Product Datasheet

COMMD1 Antibody (orb1254337)

Description	COMMD1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human COMMD1 (NP_689729.1).
Target	COMMD1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N668
Dilution Range	IF: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of L929 cell...



Immunofluorescence analysis of U2OS cell...