



Product Datasheet

LGALS3BP Antibody (orb1254356)

Description	LGALS3BP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 306-585 of human LGALS3BP (NP_005558.1).
Target	LGALS3BP
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 75kDa-110kDa
Uniprot ID	Q08380
Dilution Range	WB: 1:500 - 1:2000IF: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot
analysis of extracts
of var...



Immunofluorescence analysis of MCF7 cell...