


Product Datasheet

SETMAR Antibody (orb1254377)

Description	SETMAR Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human SETMAR (NP_006506.3).
Target	SETMAR
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q53H47
Dilution Range	WB: 1:500 - 1:2000IF: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of HeLa
cell...