



Product Datasheet

Arrestin 1 (phospho-Ser412) antibody (orb125451)

Description	Rabbit Polyclonal to Arrestin 1 (phospho-Ser412).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P49407
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from C...



Immunohistochemistry analysis of formali...