


Product Datasheet

SNRPG Antibody (orb1254528)

Description	SNRPG Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic peptide of human SNRPG
Target	SNRPG
Preservatives	PBS with 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000IF: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 10kDa
Uniprot ID	P62308
Dilution Range	WB: 1:500 - 1:1000IF: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Jur...