



Product Datasheet

m6A Antibody (orb1254705)

Description	m6A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	Nucleotide Array
Immunogen	Chemical/ Small Molecule corresponding to N6-methyladenosine (m6A) conjugated to keyhole limpet haemocyanin.
Target	m6A
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



The m6A rabbit polyclonal antibody (4 µ...



The m6A rabbit polyclonal antibody (500n...