


Product Datasheet

H4K5me1 Antibody (orb1254706)

Description	H4K5me1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	A synthetic peptide of human MonoMethyl-Histone H4-K5.
Target	H4K5me1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of NIH-3T3
C...