


Product Datasheet


MKI67 Antibody (orb1254716)

Description	MKI67 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant protein of human MKI67
Target	MKI67
Preservatives	PBS with 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46013
Dilution Range	IF: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of NIH-3T3

C...