



Product Datasheet

CASQ1 Antibody (orb1254731)


Description	CASQ1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-280 of human CASQ1 (NP_001222.3).
Target	CASQ1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IHC: 1:50 - 1:200IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31415
Dilution Range	IHC: 1:50 - 1:200IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...



Immunofluorescence analysis of NIH-3T3 c...