


Product Datasheet

CD70 Antibody (orb1254733)

Description	CD70 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant protein of human CD70.
Target	CD70
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P32970
Dilution Range	IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of rat
splee...