


Product Datasheet


ERCC5 Antibody (orb1254757)

Description	ERCC5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ERCC5 (NP_000114.2).
Target	ERCC5
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 133kDa
Uniprot ID	P28715
Dilution Range	IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of NIH-3T3

C...