


Product Datasheet

PLA2G4D Antibody (orb1254782)

Description	PLA2G4D Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-680 of human PLA2G4D (NP_828848.3).
Target	PLA2G4D
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86XP0
Dilution Range	IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of HeLa
cell...