



Product Datasheet

BrdU Antibody (orb1254806)

Description	BrdU Antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Bromodeoxyuridine
Target	BrdU
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	IHC: 1:50 - 1:500 ICC: 1:50 - 1:500 IF: 1:50 - 1:500
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC: 1:50 - 1:500 ICC: 1:50 - 1:500 IF: 1:50 - 1:500
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunofluorescence analysis of HeLa cell...