


Product Datasheet

NR4A1 Antibody (orb1254909)

Description	NR4A1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human NUR77 (NP_775180.1).
Target	NR4A1
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 75kDa
Uniprot ID	P22736
Dilution Range	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...