


Product Datasheet

SLC25A19 Antibody (orb1254922)

Description	SLC25A19 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human SLC25A19 (NP_068380.3).
Target	SLC25A19
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:2000IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HC21
Dilution Range	WB: 1:1000 - 1:2000IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of NIH/3T3
C...