


Product Datasheet

APOBEC3A Antibody (orb1254939)

Description	APOBEC3A Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human APOBEC3A (NP_663745.1).
Target	APOBEC3A
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:200 - 1:500
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 26kDa
Uniprot ID	P31941
Dilution Range	WB: 1:200 - 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HT-...