



Product Datasheet

CREB1 Antibody (orb1254959)


Description	CREB1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human CREB1 (NP_604391.1).
Target	CREB1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IHC: 1:50 - 1:200IF: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 43kDa
Uniprot ID	P16220
Dilution Range	WB: 1:500 - 1:2000IHC: 1:50 - 1:200IF: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of A-431 cells.



Immunohistochemistry of paraffin-embedded tissue sections.



Immunohistochemistry of paraffin-embedded tissue sections.