


Product Datasheet

ANXA1 Antibody (orb1255004)

Description	ANXA1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	A synthetic peptide of human ANXA1
Target	ANXA1
Preservatives	PBS with 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000IHC: 1:50 - 1:100Flow: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 34kDa
Uniprot ID	P04083
Dilution Range	WB: 1:500 - 1:1000IHC: 1:50 - 1:100Flow: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of A20...