


Product Datasheet

MST1R Antibody (orb1255020)

Description	MST1R Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human MST1R
Target	MST1R
Preservatives	PBS with 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 145kDa
Uniprot ID	Q04912
Dilution Range	WB: 1:500 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...