


Product Datasheet

CAPN12 antibody (orb125508)

Description	Rabbit Polyclonal to CAPN12.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q6ZSI9
Dilution Range	WB: 1:500-1:3000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
extracts
from N...