



Product Datasheet

CARD6 antibody (orb125510)


Description	Rabbit Polyclonal to CARD6.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9BX69
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:500-1:1000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...