


Product Datasheet

BL-CAM (phospho-Tyr807) antibody (orb125530)

Description	Rabbit Polyclonal to BL-CAM (phospho-Tyr807).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P20273
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from K...