


Product Datasheet

CDC25C (phospho-Ser198) antibody (orb125539)

Description	Rabbit Polyclonal to CDC25C (phospho-Ser198).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P30307
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...