


## Product Datasheet

### CDH26 antibody (orb125546)

|                            |                                   |
|----------------------------|-----------------------------------|
| <b>Description</b>         | Rabbit Polyclonal to CDH26.       |
| <b>Species/Host</b>        | Rabbit                            |
| <b>Reactivity</b>          | Human                             |
| <b>Conjugation</b>         | Unconjugated                      |
| <b>Tested Applications</b> | IHC                               |
| <b>Immunogen</b>           | Internal                          |
| <b>Concentration</b>       | 1mg/ml                            |
| <b>Note</b>                | For research use only             |
| <b>Clonality</b>           | Polyclonal                        |
| <b>Uniprot ID</b>          | <b>P10632</b>                     |
| <b>Dilution Range</b>      | IHC-P: 1:50-1:100, ELISA: 1:20000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Immunohistochemical analysis of formalin...