


Product Datasheet

p19 INK4d antibody (orb125556)

Description	Rabbit Polyclonal to p19 INK4d.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P55273
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...