



Product Datasheet

CALB2 Antibody (orb1255575)


Description	CALB2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-271 of human CALB2 (NP_001731.2).
Target	CALB2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:4000IHC: 1:50 - 1:100IF: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 29kDa
Uniprot ID	P22676
Dilution Range	WB: 1:1000 - 1:4000IHC: 1:50 - 1:100IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunofluorescence analysis of rat brain...