


Product Datasheet

C/EBP-epsilon (phospho-Thr74) antibody (orb125561)

Description	Rabbit Polyclonal to C/EBP-epsilon (phospho-Thr74).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q15744
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...