


## Product Datasheet

### CLIP4 antibody (orb125591)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | Rabbit Polyclonal to CLIP4.     |
| <b>Species/Host</b>        | Rabbit                          |
| <b>Reactivity</b>          | Human                           |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | WB                              |
| <b>Immunogen</b>           | Internal                        |
| <b>Concentration</b>       | 1mg/ml                          |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Polyclonal                      |
| <b>Uniprot ID</b>          | <b>Q8N3C7</b>                   |
| <b>Expiration Date</b>     | 12 months from date of receipt. |



Western blot analysis of extracts from H...