



Product Datasheet

Stefin A antibody (orb125629)


Description	Rabbit Polyclonal to Stefin A.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01040
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from A...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...