



Product Datasheet

GTF2H3 Antibody (orb1256540)


Description	GTF2H3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-308 of human GTF2H3 (NP_001507.2).
Target	GTF2H3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:50 - 1:100IP: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 36kDa
Uniprot ID	Q13889
Dilution Range	WB: 1:500 - 1:2000IF: 1:50 - 1:100IP: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of A549 cell...



Immunoprecipitation analysis of 200 ug e...