



Product Datasheet

PYGB Antibody (orb1257037)


Description	PYGB Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 544-843 of human PYGB (NP_002853.2).
Target	PYGB
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IF: 1:20 - 1:100IP: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 105kDa
Uniprot ID	P11216
Dilution Range	WB: 1:500 - 1:2000IF: 1:20 - 1:100IP: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of MCF-7 cel...



Immunoprecipitation analysis of 200 ug e...