



Product Datasheet

TNRC6A Antibody (orb1257195)

Description	TNRC6A Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 950-1050 of human TNRC6A (NP_055309.2).
Target	TNRC6A
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:1000IF: 1:50 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 180kDa
Uniprot ID	Q8NDV7
Dilution Range	WB: 1:500 - 1:1000IF: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of U2OS cell...