



## Product Datasheet

### PRPF8 Antibody (orb1257215)


|                            |                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | PRPF8 Antibody                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                 |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                           |
| <b>Tested Applications</b> | IF, IHC, IP, WB                                                                                        |
| <b>Immunogen</b>           | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PRPF8 (NP_006436.3). |
| <b>Target</b>              | PRPF8                                                                                                  |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.                                                     |
| <b>Form/Appearance</b>     | Liquid                                                                                                 |
| <b>Concentration</b>       | batch dependent                                                                                        |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                                                            |
| <b>Note</b>                | For research use only                                                                                  |
| <b>Application notes</b>   | WB: 1:500 - 1:1000IHC: 1:50 - 1:100IF: 1:50 - 1:100IP: 1:50 - 1:100                                    |
| <b>Isotype</b>             | IgG                                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                             |
| <b>MW</b>                  | Observed: 274kDa                                                                                       |
| <b>Uniprot ID</b>          | <b>Q6P2Q9</b>                                                                                          |
| <b>Dilution Range</b>      | WB: 1:500 - 1:1000IHC: 1:50 - 1:100IF: 1:50 - 1:100IP: 1:50 - 1:100                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                        |



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...