


Product Datasheet

FAM3D antibody (orb125730)

Description	Rabbit Polyclonal to FAM3D.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96BQ1
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
extracts
from K...