



Product Datasheet

TP53 Antibody (orb1257343)


Description	TP53 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IHC, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 218-390 of mouse p53 (NP_035770.2).
Target	TP53
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:3000 - 1:6000IHC: 1:50 - 1:200IP: 1:20 - 1:50ChIP: 1:20 - 1:50
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 53kDa
Dilution Range	WB: 1:3000 - 1:6000IHC: 1:50 - 1:200IP: 1:20 - 1:50ChIP: 1:20 - 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...