



Product Datasheet

PRDM14 Antibody (orb1257447)


Description	PRDM14 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-571 of human PRDM14 (AAH52311.1).
Target	PRDM14
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IHC: 1:50 - 1:200IF: 1:50 - 1:200IP: 1:20 - 1:100ChIP: 1:20 - 1:100
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 64kDa
Uniprot ID	Q9GZV8
Dilution Range	WB: 1:500 - 1:2000IHC: 1:50 - 1:200IF: 1:50 - 1:200IP: 1:20 - 1:100ChIP: 1:20 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of rat...



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of NIH/3T3 C...